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Suction pump massage – SPM – is a combination of cupping and massage (cupping massage). It is a holistically oriented, cosmetic treatment for the face or whole body. It is also suitable for applications in physiotherapy.

SPM

Suction pump massage

In SPM a specific pre-selected negative pressure is conducted via connecting tubes into treatment pistols which are then used to massage the customer. **Different pistols are available for different areas of the body.** It is possible to work with two pistols at the same time.

Depending on the intensity of the negative pressure and the massage technique, a targeted lymph drainage can be carried out with SPM very simply, skin tissue regenerated, fat metabolism stimulated (ideal for fatty deposits or jodhpur thighs) or the tissue firmed. A further speciality is the possibility of phagocytosis.

Suction pump massage usually leads to **immediately visible and noticeable results.** It is suitable for all customers who would like a fast, long-lasting result. With regard to facial and décolleté care, SPM is suitable not only for regeneration (anti-ageing), lifting and wrinkle treatment but also for acne and rosacea. In the whole body area SPM meets the requirements for **Differentiated Cellulite Treatment** in terms of lymph drainage, activation of the fat metabolism and connective tissue firming. SPM also provides the possibility of special treatment for pregnancy marks. It is also ideal for relaxing massages or decongesting cramped muscles.


Treatments with SPM take between 30 and 60 minutes depending on the treatment aim and the treatment zone. In cosmetic facial care SPM is also often **incorporated as a component of existing treatment plans** (e.g. as a preparatory lymph drainage, for regeneration or as lifting). A particular advantage is the possibility of working with two pistols at the same time which reduces the treatment times compared to working with single applicators. The treatment pressure is adjustable between 0 and 500 mbar (-). Option of permanent or pulsating negative pressure. The pistols are made of special glass, are easy to clean and always guarantee optimum hygiene.

Technical details

Mains voltage	230 V / 50-60 Hz	115 V / 50-60 Hz
Current intake	490 mA	1000 mA
Weight (device)	7.5 kg	7.5 kg
Weight (accessories)	3.5 kg	3.5 kg
Size WDH (device)	450 x 400 x 190 mm	450 x 400 x 190 mm

Between 40 and 90 euro is paid for SPM treatments depending on the duration and zone. About 15 euro is added to the basic charge for use as a component in facial care.

Delivered as a **complete package** consisting of control unit, eight pistols in four different sizes, connecting tubes, detailed user manual, overview poster for the application technique in the cubicle and advertising material.

- Certified quality according to ISO / EN 9001
- Made in Germany
- 2-year guarantee
- Guarantee extension possible
- Also available as a medical device according to ISO 13485

Ref: 24202