



Dermio Care[®] Plus

Plasma and color light for problem skin,
anti aging and revitalization


Innovative

„Dermio Care® Plus“* is an innovative method for cosmetic and medical treatments. During the application the customers are treated with the energy of oxygen in the form of applied electrically negatively charged ions (Plasma treatment). The application pursues several goals: visible skin rejuvenation already after the first use, treatment of problem skin, maintenance of skin health, increase of vitality and positive support of the metabolism in the entire organism.

The Dermio Care® Plus treatment is based on findings of quantum physics and is patented. Negative ions (plasma) are gained from the ambient air. The energy is directed into the specially coated „Top Head“ where it is combined with the color light therapy. The ions penetrate through the skin to promote the cell metabolism. At the same time the energy is inhaled through the mouth and nose and is immediately available for the whole organism, all organs and the muscles.

* dermis = skin / ion = electrical particle, "plus" stands for the additional phototherapy module.

The more the better. Our organism needs the electricity of the negative ions from the oxygen molecule.



Positive and negative ions have different influences on our organism.

Affected Target	Blood	Heart	Blood Pressure	Respiratory System	Blood Vessels	Bones	Vitality
When rich in negatively ions	Weak Alkaline pH balance	Works efficiently	Normal	Smooth	Expanded	Healthy	Vigorous
When rich in positive ions	Acidic pH Balance	Has to work hard	High	Difficult	Constricted	Weak	Fatigued

The Dermio Care® Plus application is aimed at all customers with high expectations in anti-aging treatments and immediately visible and noticeable results.

Color light therapy: scientifically proven and effective

Color light treatment (phototherapy) has been known for a long time and has been used successfully in medicine for decades. There was even a Nobel Prize **. The effect of the colored light is reinforced by the negative ions.



Applications of red light therapy:

- Wound healing
- Increase in collagen formation / collagen synthesis ***
- Skin repair (anti-aging)
- Prevents and reduces scar tissue
- Stimulates ATP in cells and increases cellular reproduction



Applications of green light therapy:

- Detoxification
- Antibacterial
- Relieves stress
- Strengthens the immune system
- Builds muscle and bones
- As aphrodisiac and sexual tonic



Applications of blue light therapy:

- Eczema / psoriasis
- Acne
- Antiseptic and antibacterial
- Anti-inflammatory and astringent effect
- Anti sebum production
- Soothes the nerves
- Strengthens the immune system

** Nobel Prize for Medicine, Niels Ryberg Finsen, 1903

*** See study by Wish, A.: A Controlled Trial to Determine the Efficacy of Red and Near-Infrared Light Treatment in Patient Satisfaction, Reduction of Fine Lines, Wrinkles, Skin Roughness, and Intradermal Collagen Density Increase. In: Photomedicine and Laser Surgery, 2014 (2), 93-100

Effective and profitable

The innovative application is carried out conveniently while sitting or lying down. The regular treatment duration is 30 minutes.

In cosmetic facial care Dermio Care® Plus is suitable both before or after regular care as well as within the appointment during mask break (also in combination with preparations, e.g. the Crystal Collagen Gel Mask from Weyergans high care® cosmetics). In this case, the treatment duration is reduced by 15 minutes.

Dermio Care® Plus is also optimally suitable in connection with a full-body application (e.g. gliding stroke massage, body wraps, Vacustyler®, ultrasound or dry water massage).

The Dermio Care® Plus applications can also be sold as a cure, e.g. as a 10x cure distributed over a skin regeneration cycle (four weeks). Also ideal as a half year or annual cure or as a module in connection with a membership.

Additional advantages:

- ✓ Ideal for every holistic cosmetic institute
- ✓ Easy to apply (application treatment)
- ✓ No staffing
- ✓ No wear and tear, no subsequent costs
- ✓ Suitable for everyone, no contraindications
- ✓ Low costs (from 137 Euro per month)
- ✓ Complete marketing and advertising package
- ✓ Single treatment price 39,80 Euro / 10x card only 280 Euro
- ✓ CE mark (TÜV Hessen) and Made in Germany
- ✓ Future-proof
- ✓ Two-year warranty (renewable)

Technical specifications

Voltage	230 V / 45 VA
Weight	9,0 kg
Dimensions (W D H)	450 x 420 x 140 mm